

Fractionation Of Milk For Trace Analysis Of Contaminants And Residues

Christian BAUMGARTNER¹ and Jan-Michael STEILS¹

Aim

Developing a new **passive method** for **raw milk** to prepare samples for **trace analyses of residues** of chlorate and perchlorate (C/PC) with **IC** and **LC-MS/MS**.

Materials and Methods

Unprocessed raw milk was separated into die phases „water" and „fat & protein,, with a fractionation unit (FraMiTrACR C/PC[®]) and a standard benchtop centrifuge. Both, 10 mL raw milk using a fixed angle rotor at 30 minutes and 2,000 x g as well as 5 mL raw milk using a "swing out" rotor at 30 minutes and 4,000 x g resulted in half of the volume of the feed added recovered as the "water phase".

For the analysis of the water phase, the **930 Compact IC Flex system (Metrohm)** with Dosis gradient technique was used for the determination of anions after sequential suppression and conductivity detection. For each determination, **0.25 mL of water phase** was injected into the analyser. At the beginning of the development, the anion matrix in the water phase was collected from milk samples that were free of chlorate and perchlorate. To improve the detection limit, the characteristic anion matrix of the water phase has been subtracted when evaluating the results.

The chlorate and perchlorate content was determined using **spiked samples in standard series** and also in **native samples**. In order to obtain **comparability** with previously used methods, the investigated raw milk samples were simultaneously analysed in a **contract laboratory** using the **modified Quppe method** and liquid chromatography with tandem mass spectrometry coupling (LC-MS/MS).

Results

In the **water phase**, prepared with the fractionation unit (FraMiTrACR C/PC[®]) a **detection limit of 0.003 mg/kg chlorate and perchlorate** could be achieved by using the described **IC method**. Using a **LC-MS/MS method**, a **detection limit of down to 0.001 mg/kg chlorate and perchlorate** could be achieved.



In direct **comparison** to well-established sample preparation methods, such as **Quechers** and **Quppe**, there is a **significant saving in working time (about 70%)** and thus in personell costs as well.

In addition, the complete **elimination of extraction solvents** or other additives leads to further cost savings in the laboratory process chain, simultaneously **excluding the risk of unwanted contamination**.

The **fractionation unit FraMiTrACR[®]** offers in various specifications **new opportunities** for trace analysis in milk and milk products.

Positive results for **antibiotics** are available as well.

Further studies for other analytes will be following soon.



References

1. Chlorate and Other Oxychlorine Contaminants Within the Dairy Supply Chain, 2018 Institute of Food Technologists[®], doi: 10.1111/1541-4337.12393, Vol.17,2018 Comprehensive Reviews in Food Science and Food Safety, William P. McCarthy, Tom F. O'Callaghan, Martin Danahar, David Gleeson, Christine O'Connor, Mark A. Fenelon and John T. Tobin
2. Perchlorate in Milk, Environ. Sci. Technol. 2003, 37, 21, 4979–4981, Andrea B. Kirk, Ernest E. Smith, Kang Tian, Todd A. Anderson, and Purnendu K. Dasgupta
3. Quick Method for the Analysis of Numerous Highly Polar Pesticides in Food Involving Extraction with Acidified Methanol and LC MS/MS Measurement II. Food of Animal Origin (QuPpe AO Method), M. Anastassiades; A.-K. Wachtler; D. I. Kolberg; E. Eichhorn; A. Benkenstein; S. Zechmann; D. Mack; A. Barth; C. Wildgrube; I. Sigalov; S. Görlich; D. Dörk; G. Cerchia
4. Sample processing method for the determination of perchlorate in milk, Analytica Chimica Acta 567 (2006) 73–78, Jason V. Dyke, Andrea B. Kirk, P. Kalyani Martinelango, Purnendu K. Dasgupta
5. DIN ISO 11352:2013 Wasserbeschaffenheit - Abschätzung der Messunsicherheit beruhend auf Validierungs- und Kontrolldaten (ISO 11352:2012)



¹ pureMilk analytical GmbH

www.pureMilk.eu

Address for correspondence:

Dr. Christian Baumgartner

Freyunger Str. 15 | D-94065 Waldkirchen

office@pureMilk.eu | fon +49 171 8462612